

**INNOVACYN, INC.**  
**3546 N. RIVERSIDE AVE., RIALTO, CA, 92377**

Toll-Free: 866-318-3116

Fax: 909-428-1947

Website: [www.innovacyn.com](http://www.innovacyn.com)

Website: [www.vetericyn.com](http://www.vetericyn.com)

Email: [customerservice@innovacyn.com](mailto:customerservice@innovacyn.com)

---

Every effort has been made to ensure the accuracy of the information published. However, it remains the responsibility of the readers to familiarize themselves with the product information contained on the USA product label or package insert.

---

## **VETERICYN PLUS<sup>®</sup> VF ANTIMICROBIAL WOUND & SKIN HYDROGEL**

**Innovacyn**

**SPRAYABLE HYDROGEL**

Wound Management / Moisture Donating

- Acute & Chronic Wounds
- Post-Surgical Applications
- Abrasions
- Burns

**HYPOCHLOROUS TECHNOLOGY**

**non-sensitizing**

**non-cytotoxic**

**non-irritating**

**Dispensed by licensed veterinarians only.**

**Indications:** Intended for use by veterinary professionals for moistening and debriding wounds of all types, burns, pressure sores, ulcers, abrasions, dermal and mucosal irritations, post-surgical incisions, grafted and donor sites, in addition to moistening and lubricating absorbent wound dressings. Vetericyn Plus<sup>®</sup> Antimicrobial Wound & Skin Hydrogel relieves itch and pain often associated with dermal irritations. Promotes a moist environment which has been proven to help wounds heal faster.

**Directions for use: Shake well before each use.** Remove excess hair and necrotic tissue from wound area. Necrotic tissue can be softened with hydrogel for easier debridement. Remove cap and saturate the affected area. Reapply as needed to keep wound moist but not wet. If wound requires a bandage, apply at each bandage change. For best results use Vetericyn Plus<sup>®</sup> Wound & Skin Cleanser to flush and clean the area prior to applying hydrogel. Faster healing is achieved when wounds are clean and maintained in a moist environment. Safe for use around eyes, ears, nose and mouth.

**Use Precautions:** For external use only, not for injection. For animal use only, not for human use. Keep out of reach of children. Always consult with a veterinarian for systemic conditions. Discontinue use if irritation develops and contact your veterinarian.

**Storage:** Store at room temperature away from direct sunlight and heat. Do not allow to freeze. Do not use after expiration date. After use, seal or close dispenser.

**Stability:** See Expiration Date

**Disposal:** No special disposal requirements. Vetericyn<sup>®</sup> Plus VF is environmentally friendly.

**Ingredients: Active:** Hypochlorous Acid (0.012%) **Inactive:** Electrolyzed Water, Sodium Chloride, Sodium Hypochlorite, Lithium Magnesium Sodium Silicate, Sodium Bicarbonate, Phosphates, Sodium Sulfate

Innovacyn products are manufactured under strict quality guidelines of ISO 13485 and FDA QSR. Vetericyn Plus hydrogel formula contains Hypochlorous Acid, a known antimicrobial, which serves as a preservative to inhibit the growth of microorganisms in the solution.

Manufactured & Bottled by Innovacyn, Inc., 3546 N. Riverside Ave, Rialto CA 92377

**Made in the USA**

**vetericynVF.com**

**866.318.3116**

	Part Number:
2 fl oz	2042; 20042 v3.0
16.9 fl oz (500 ml)	2040; 90021 v5.0

**CPN:** 1433006.8

---

By using this content, you agree to the [Terms of Service](#) and [Privacy Policy](#).

Copyright 2020 - Animalytix LLC